

INSTRUCTIONS FOR USE



COMFORTEAR® LACRISOLVE™ 180 ABSORBABLE PUNCTUM PLUGS FOR MEDIUM-TERM OCCLUSION THERAPY (180 DAYS)

INTRODUCTION

Comfortear® Lacrisolve™ 180 Absorbable Punctum Plugs provide medium-term occlusion therapy of the tear drainage system lasting approximately 180 days.

- Are effective in the horizontal canaliculi
- Never touch the eye
- Are comfortable after proper placement, and will not fall out of the punctum

Punctal dilation and topical anesthesia are not typically required for the insertion of the Comfortear® Lacrisolve™ 180 Absorbable Punctum Plugs.

INDICATIONS

The Comfortear® Lacrisolve™ 180 Absorbable Punctum Plugs are intended to temporarily block tear drainage by the occlusion of the canaliculus in order

to:

- Determine the potential effectiveness of permanent occlusion
- Temporarily treat dry eye syndrome, and the dry eye components of various ocular surface diseases
- Temporarily enhance the efficacy of topical medications or ocular lubricants
- Temporarily treat contact lens intolerance secondary to dry eye
- Temporarily treat dry eye after ocular surgery

CONTRAINDICATIONS

Comfortear® Lacrisolve™ 180 Absorbable Punctum Plugs are contraindicated for use in patients with known sensitivity to any one of the materials that comprise the Comfortear® Lacrisolve™ 180 Absorbable Punctum Plugs, infectious conjunctivitis, dacryocystitis, inflammation of the eyelid, infected eyes, or epiphora.

BEFORE INSERTION

Patients should be given a complete eye examination with thorough history. Patients with pre-existing intermittent tearing should receive pressure irrigation to rule out canalicular obstruction.

PACKAGING

Comfortear® Lacrisolve™ 180 Absorbable Punctum Plugs for medium-term occlusion therapy are made of Polydioxanone and provided sterile (by EO Sterilization), two plugs per pouch, one pouch per box.

Plugs are available in two sizes: 0.4mm and 0.5 mm.
Insert the largest size possible for maximum occlusion benefits.

This device is designed, intended and distributed for single use only. DO NOT RE-STERILIZE OR REUSE THIS DEVICE. There is no data to support sterility and functionality of the device after reprocessing.

MADE IN USA (CE 0470)

CP Medical, Inc. 1775 Corporate Dr, Ste 150, Norcross, GA 30093 USA



www.altacor-pharma.com

For further details please call +44(0)118 902 6766